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

1. Markwick A, Lee L, Horsfall M, Sinhal A, Chew D. TCT-784 Prognostic implications of moderate and severe mitral regurgitation in contemporary clinical care [abstract]. J Am Coll Cardiol. 2012;60(17)(suppl B):B228. 2. Grigioni F, Tribouilloy C, Avierinos JF, et al; MIDA Investigators. Outcomes in mitral regurgitation due to flail leaflets: a multicenter European study. JACC Cardiovasc Imaging. 2008;1(2):133-141. 3. Enriquez-Sarano M, Avierinos JF, Messika-Zeitoun D, et al. Quantitative determinants of the outcome of asymptomatic mitral regurgitation. N Engl J Med. 2005;352(9):875-883. 4. Bursi F, Barbieri A, Grigioni F. Prognostic implications of functional mitral regurgitation according to the severity of the underlying chronic heart failure: a long-term outcome study. Eur J Heart Fail. 2010;12(4):382-388. 5. Ahmed A, Aronow WS, Fleg JL. Higher New York Heart Association classes and increased mortality and hospitalization in patients with heart failure and preserved left ventricular function. Am Heart J. 2006;151(2): 444-450. 6. Montant P, Chenot F, Robert A, et al. Long-term survival in asymptomatic patients with severe degenerative mitral regurgitation: a propensity score-based comparison between an early surgical strategy and a conservative treatment approach. J Thorac Cardiovasc Surg. 2009;138(6):1339-1348.

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1. Nkomo VT, Gardin JM, Skelton TN, Gottdiener JS, Scott CG, Enriquez-Sarano M. Burden of valvular heart diseases: a population-based study. *Lancet*. 2006;368(9540):1005-1011.
2. United States Census Bureau, 2010. Age and sex composition. Issued May 2011:1-16.
3. Abbott data on file as of March 2019.

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
	Through 1 Year # (%) of patients	1 Year to 5 Years # (%) of patients
Single Leaflet Device Attachment (SLDA)	10 (6.3%)	0 (0.0%)
MV stenosis	1 (0.6%)	0 (0.0%)
Device Embolization	0 (0.0%)	0 (0.0%)


Based on N=158 who were implanted with 1 or 2 MitraClip devices

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

1. Feldman T, Kar S, Elmariah S, et al. Randomized comparison of percutaneous repair and surgery for mitral regurgitation: 5-year results of EVEREST II. J Am Coll Cardiol. 2015;66(25):2844-2854.

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COAPT™ TRIAL

Stone GW, Lindenfeld J, Abraham WT, et al.

Transcatheter mitral-valve repair in patients with heart failure. N Engl J Med. September

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Recommendations for Intervention for Chronic Primary MR/器質性MR

ACC/AHA CLINICAL PRACTICE GUIDELINE


2020 ACC/AHA Guideline for the Management of Patients With Valvular Heart Disease

A Report of the American College of Cardiology/American Heart Association Joint Committee on Clinical Practice Guidelines

推奨クラス	エビデンスレベル	
2a	B-NR	6. In severely symptomatic patients (NYHA class III or IV) with primary severe MR and high or prohibitive surgical risk, transcatheter edge-to-edge repair (TEER) is reasonable if mitral valve anatomy is favorable for the repair procedure and patient life expectancy is at least 1 year. ^{17,18}

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Recommendations for Intervention for Chronic Secondary MR/機能性MR

ACC/AHA CLINICAL PRACTICE GUIDELINE

2020 ACC/AHA Guideline for the Management of Patients With Valvular Heart Disease


A Report of the American College of Cardiology/American Heart Association Joint Committee on Clinical Practice Guidelines

推奨クラス	エビデンスレベル	
2a	B-R	1. In patients with chronic severe secondary MR related to LV systolic dysfunction (LVEF <50%) who have persistent symptoms (NYHA class II, III, or IV) while on optimal GDMT for HF (Stage D), TEER is reasonable in patients with appropriate anatomy as defined on TEE and with LVEF between 20% and 50%, LVESD ≤70 mm, and pulmonary artery systolic pressure ≤70 mm Hg. ¹⁻⁸

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JCS/JATS/JSVS/JSCS 2020 Guideline on the Management of Valvular Heart Disease

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MR may develop in patients with preserved LV systolic function who have progressive LA dilation, leading to enlargement of the mitral annulus and malcoaptation of the leaflets. This may arise in conditions such as HF with preserved LVEF, restrictive cardiomyopathy, and nonobstructive hypertrophic cardiomyopathy. These patients often have associated AF, which may contribute to the progression of LA and annular dilation, thus increasing the severity of MR and successful ablation of AF may reduce or eliminate MR.

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MitraClip is a trademark of the Abbott Group of Companies.

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Primary Efficacy Endpoint
Acute Procedural Success

Implant Rate.....100% (30/30)

Acute Procedural Success86.7% (26/30)

- Four subjects were non-APS due to MR>2+ at discharge
 - MR Severity for non-APS patients remained the same or improved
 - NYHA Classification for non-APS patients remained the same or improved
 - Procedure and clip did not appear to exacerbate the mitral insufficiency or clinical condition

Primary Safety Endpoint
MAE at 30 Days and 1 Year

Category of MAE	Through 30 Days (CEC-Adjudicated)		31 Days Through 1 Year (Site-Reported)	
	% Patients (n/N)	# Events	% Patients (n/N)	# Events
Death	0.0 (0/30)	0	3.3% (1/30)	1
Stroke	0.0 (0/30)	0	0.0 (0/30)	0
Myocardial Infarction	0.0 (0/30)	0	0.0 (0/30)	0
Renal Failure	0.0 (0/30)	0	3.3% (1/30)	1
Non-elective cardiovascular surgery for device or procedure related adverse events	0.0 (0/30)	0	0.0 (0/30)	0
MAE rate (%)	0.0 (0/30)	0	6.7% (2/30)	2

Note: All data through 44 and 393 days post index procedure are included for the 30-Day and 1-Year MAE rate, respectively

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

*Voluntary recall of MitraClip for 6 months in 2011.

Data source: German Federal Statistics Office databases (DRG statistics) for TMVr and mitral valve surgery for the period 2010 through 2014.

Proprietary and confidential — do not distribute

1. Pietzsch JB, Weber SA, Pietzsch ML, Conradi L. Adoption of transcatheter mitral valve repair in Germany: Utilization patterns and case volumes compared to mitral valve surgery in the period 2010-2014. Value in Health. 2016;19:A356.

2. Downs E, Lim S, Ragosta M, et al. The impact of a percutaneous mitral repair program on surgical mitral valve volume. J Thorac Cardiovasc Surg. 2015;150(5):1093-1097.

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STRUCTURAL INTERVENTIONS

Occlusion therapies, including congenital heart and stroke-risk reduction

PATENT DUCTUS ARTERIOSUS

**ATRIAL SEPTAL DEFECTS AND
PATENT FORAMEN OVALE**

**VENTRICULAR
SEPTAL DEFECTS**

VALVULAR THERAPIES


Transcatheter and open-surgical valve repair and replacement

MITRAL VALVE DISEASE

AORTIC VALVE DISEASE

AMJ STRUCTURAL HEART 取扱製品

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40 YEARS

SURGICAL MARKET LEADERSHIP

Mechanical Valves

>2 million implanted

Tissue Valves

MASTERS™

EPIC™ 1.0

16+ YEARS

REAL-WORLD TRANSCATHETER EXPERIENCE

Transcatheter Mitral Valve Repair

100K patients treated

MITRACLIP™

MITRACLIP™ NT
(2016)

MITRACLIP™ NTR/XTR
(2018)

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MITRACLIP™ G4+
(2019)

AMJ STRUCTURAL HEART 取扱製品

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40 YEARS

SURGICAL MARKET LEADERSHIP

Mechanical Valves
>2 million implanted

Tissue Valves

REGENT™


EPIC™ SUPRA

TRIFECTA™ GT

AMJ STRUCTURAL HEART 取扱製品

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20 YEARS

LEADERSHIP BEYOND QUESTION FOR CONFIDENCE BEYOND MEASURE

*Atrial Septal
Defects
AMPLATZERTM
SEPTAL
OCCLUDERS*

*Patent
Ductus
Arteriosus
AMPLATZERTM
DUCT
OCCLUDERS*

*Patent
Foramen Ovale
AMPLATZERTM
PFO OCCLUDER*

*Patent Ductus
Arteriosus
AMPLATZERTM
PICCOLOTM
OCCLUDER*

*Delivery System
AMPLATZERTM
TREVISIOTM
SYSTEM**

Proprietary and confidential — do not distribute

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AMJ STRUCTURAL HEART 取扱製品

PROVIDING PATIENTS LEADING THERAPIES

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©2021 Abbott. All rights reserved. (APJ00001202-JPN-Rev.A)